



## 2007 Battery Life Expectancy

Standard Products	Life Expectancy (in days)	Batteries Required (in 12 months)
<b>SIZE 675</b>		
DaVinci P X P	22	17
<b>SIZE 13</b>		
Destiny BTE	11	33
Axent BTE	13	28
S13 Arista BTE	10	37
J13 Sequel BTE	16 - 20	18 - 23
J13 Endeavour BTE	16 - 20	18 - 23
J13 Mesa	16	23
<b>SIZE 312</b>		
AristaLite BTE	12.5	30
<b>SIZE 10</b>		
Aspect OTE	8	46

Custom Products	Life Expectancy (in days)	Batteries Required (in 12 months)
<b>SIZE 13</b>		
Destiny ITE	12 - 16	23 - 31
Axent Custom	23	16
Arista Custom	23	16
Sequel Custom	22	17
Endeavour Custom	22	17
Mesa	19 - 21	17 - 19
Cierra	19 - 21	17 - 19
Discovery D	16 - 24	15 - 23
Inteli D	23 - 24	15 - 16
Inteli Power D	15 - 23	16 - 25
<b>SIZE 312</b>		
Destiny ITC/ITE	6 - 9	41 - 61
Axent Custom	8 - 12	31 - 46
Arista Custom	8 - 12	31 - 46
Sequel Custom	8 - 11	33 - 46
Endeavour Custom	8 - 11	33 - 46
Mesa	10 - 11	33 - 37
Cierra	7 - 11	33 - 52
Discovery D	10 - 14	26 - 37
Inteli D	10 - 14	26 - 37
Inteli Power D	9 - 13	28 - 41
<b>SIZE 10</b>		
Destiny ITC/CIC	3 - 5	73 - 122
Axent Custom	5 - 7	52 - 73
Arista Custom	5 - 7	52 - 73
Sequel Custom	5 - 7	52 - 73
Endeavor Custom	5 - 7	52 - 73
Mesa	5 - 6	61 - 73
Cierra	5 - 6	61 - 73
Discovery D	5 - 8	46 - 73
Inteli D	6 - 8	46 - 61
Inteli Power D	5 - 7	51 - 73